

Ring Fluorescent Lamp for Stereo Microscope

Our Fluorescent Lamp Illuminator offers shadow-free and non-Directional illumination by fitting lamp onto microscope objective lenses. It can be used with all kinds of Microscopes. It is ideal lighting for long working without eye fatigue. Color temperatures, RED / GREEN / BLUE / Yellow filter glass made, YEX / Blacklight are available for special requirements. Application (for industrial purpose) IC Inspection, Jewel Inspection, and Production control by ITV, (for science) physiological / pathological inspection agricultural research




勝特力材料 886-3-5753170
勝特力电子(上海) 86-21-54151736
勝特力电子(深圳) 86-755-83298787
[Http://www.100y.com.tw](http://www.100y.com.tw)


The type of standard lamp and its characteristics

Lamp Type No.	Color Temperature K	Spectrum Distribution A	Peak Spectrum A	Application
80B-WW 92B-WW WarmWhite	3200 3950	4000 7200	5800	Standard White for usual application
80B-DL 92B-DL DayLight	5300 8000	4000 7200	4600	Bluish color, fit to be used for inspection of metal surface

Warm White (WW)



Daylight (DL)




80B-NEX, WW, DL 及び 各カラー色

80B-NEX, WW, DL and other colors

訂正記号	日付	訂正理由	訂正箇所	担当

勝特力材料 886-3-5753170
 胜特力电子(上海) 86-21-54151736
 胜特力电子(深圳) 86-755-83298787
[Http://www.100y.com.tw](http://www.100y.com.tw)



ピン位置、寸法はJIS.C7709.G10

編号	設計	製造	材質	カラー	重量	型式	図番
澤田	藤井	岡田	処理		個数 1	80B-NEX, WW DL, 各カラー色	RL080B00
99/7	88/6	020625	—◆—	尺度	1/1	品名	電通産業株式会社



RAYTRONICS CORP.
電通産業株式会社

勝特力材料 886-3-5753170
勝特力电子(上海) 86-21-54151736
勝特力电子(深圳) 86-755-83298787
[Http://www.100y.com.tw](http://www.100y.com.tw)




NEX



BL



YEX



REX



GEX



BEX