
100Y.COM.TW **PH-mV-TEMP. METER** Model: PH-206 *ISO-9001, CE, IEC1010*



勝特力材料 886-3-5753170 胜特力电子(上海) 86-21-34970699 WWW.100Y.COM.TW 胜特力电子(深圳) 86-755-83298787 Http://www. 100y. com. tw





The Art of Measurement

OM.TW

DIGITAL PH/mV/TEMP METER Model: PH-206 100Y.COM WEN

GENERAL SPECIFICATIONS		
Features	 * Multi functions: PH, mV, Temp. * Automatic/manual temp. compensation. * Large size (18 mm) LCD display. * Built in analog output terminal. * mV function for ion-selective, ORP and other precise mV measurement. 	
Display	18 mm (0.7") LCD, 3 1/2 digits.	
Measurement	PH, mV, Temp. (optional probe).	
Range/	PH : 0 to 14 pH x 0.01 pH	
Resolution	Temperature : 0 to 100 °C x 0.1 °C.	
Accuracy	PH : ± (0.01 pH + 1d) TEMP. : ± 0.8 ℃ mV : ± (0.5 % + 1d) <i>Temp., mV - within 23± 5 ℃ PH : Main instrument only.</i>	
Input Impedance	10^12 ohms.	
PH range Temp. Compensation	Manual : 0 $^{\circ}$ C to 100 $^{\circ}$ C external knob, adjust temp. values showed on display Auto temp. compensation : with the optional temp. probe, 0 to 100 $^{\circ}$ C.	
Analog Output	The analog output is 0.1 mV/per digit of the display both of PH, mV, Temp. function.	
Sampling Time	Approx. 0.4 second.	
Operating Temp.	0℃ to 50℃ (32°F to 122°F)	
Operating Humidity	Less than 80% RH.	
Calibration Knob	External PH 4 (slope adj.) & PH 7 (Cal. adj.).	
Power Supply	006P DC 9V battery.	
Power Consumption	Approx. 2.7 mA.	
Dimension	160 x 80 x 35 mm, (6.3 x 3.2 x 1.4 inch).	
Weight	250 g/0.55 LB (including battery).	
PH electrode (optional)	Any connector PH Electrode with BNC connector.	
Standard Accessories	Instruction Manual1 PC.	
Optional Accessories	PH electrode, ATC Temp. probe,	

APPLICATIONS					
* AQUARIUMS	* PAPER INDUSTRY	* FOOD PROCESSING	* SCHOOLS & COLLEGES		
* BEVERAGE	* PLATING INDUSTRY	* PHOTOGRAPHY	* SWIMMING POOLS		
* FISH HATCHERIES	* QUALITY CONTROL	* LABORATORY	* WATER CONDITIONING		

	OPTIONAL ACCESSORIES	1002.0	
PH ELECTRODE	General purpose pH electrode with BNC conector.		
PE-03	Range : 1 to 13 pH (typical 0 to 14 pH).	CONT. MAN MAN TO ON	
	Expoxy body, Body size : 12 mm dia. x 160 mm.	COM'L	
WV	Cable length : 1 meter.	M.TW W.	
PH ELECTRODE	General purpose & high quality pH electrode with BNC conector.		
PE-11	Range : 1 to 13 pH (typical 0 to 14 pH).		
	Expoxy body, Body size : 9.5 mm dia. x 120 mm.		
	Cable length : 1 meter.		
SPEAR PH Electrode	The "Spear Tip pH electrode " is perfect for those pH measurements in		
PH-06HD, PH-04HD	applications where sample piercing is required. Meat, sausage and cheese are		
	ideal applications. The electrode features a very durable glass measuring		
	spear packaged in a rugged virtually unbreakable epoxy body.		
	Range : 0 to 14 pH (PE-04HD), Range : 1 to 13 pH (PE-06HD)		
BUFFER SOLUTION, PH-07	PH 7.00 standard buffer solution.		
BUFFER SOLUTION, PH-04	PH 4.00 standard buffer solution.		
TEMP. PROBE, TP-06	The optional temperature probe for PH-206.	勝特力材料 886-3-5753170	
CARRYING CASE, CA-03	Vinyl soft case for PH-206. 胜特力电子(上海) 86-21-34		
CARRYING CASE, CA-06	Hard carrying case for PH-206.	胜特力电子(深圳) 86-755-8329878	

* Appearance and specifications listed in this brochure are subject to change without notice.

699 8787 Http://www. 100y. com. tw